

# MagPro R20

The optimal choice when running limited number of rTMS protocols per day

## Features and Benefits

- Low cost basic rTMS stimulator
- Highly mobile solution
- Can run most protocols for research within psychiatry
- Quick and simple protocol set-up on stimulator
- R20 Configurator for easy protocol set-up on external computer
  - Wireless/USB upload on R20
  - Possible to "lock" the R20 to a pre-defined protocol
- Optional reporting tool



## Technical specifications

Maximum Repetition Rate: 20 pps  
 Pulse mode: Standard  
 Waveform: Biphasic - 280

Output power roll-off		Mechanical data	
Max. output at 5 pps	100 %	HxWxD	150 x 390 x 440mm
Max. output at 10 pps	75 %	Weight	20kg
Max. output at 15 pps	40 %		
Max. output at 20 pps	35 %		

The choice of repetition rate and number of pulses in train puts a limit on the output power and defines a minimum Inter Train Interval.

## Compatible coils

- Static cooled coils: MCF-B65, MCF-125, MCF-B70 (up to 5500 pulses before warm-up\*)
- Non-cooled coils: C-100, C-B60, MMC-140-II, RT-120-II (up to 1500 pulses before warm-up\*)

\*Depending on the coil